

ENGG6400 Mobile Devices App Development



Group 5 Zimo Zhou (1189792), Yaowen Mei (1177855)



What is OneScan

- **OneScan** is an App is to help an auditor to conduct a more efficient and more convenient audit for a medical device factory or automotive parts factory.
- The auditor will need to open this app, scan an auditee's equipment with back camera, then with machine learning technology, the app will show the auditor
- the name of this equipment,
- the standard operation procedure of this equipment, and
- historical calibration record of this equipment.





Why OneScan

Manufactories in North America are under strict regulation:

- Automotive factories need to be audited by BIQS (IATF 16949-ISO 9001) for every 6 months if they make parts for GM and Ford,
- Medical device factories need to be inspected at least every 2 years by a FDA inspector (ISO-13485) if they make medical devices.

Calibration is one of the auditor's most interest questions[1]. The auditor's time is very precious while the auditor's knowledge is limited. Without knowing much details in the auditee's Lab, our app (Machine learning + Database) can easily provide the auditor a lot of useful information:

- 1. What is the name of this device,
- 2. Has this device been calibrated before use,
- 3. How often should it be calibrated according to the standard operation procedure,
- 4. Where can I find the historical calibration record.


Reference:


[1] Services, R. S. C. (2018, July 6). *Guide to FDA Requirements and Importance of Medical Device Calibration*. Tech Briefs. https://www.medicaldesignbriefs.com/component/content/article/mdb/tech-briefs/29754.



Mock-up --- Home Screen


- 1. User can take a picture with this app
- 2. Or choose a photo from the album
- 3. Or see the instructions about how the app is going to work





Mock-up --- 1) Take a picture

- If the device is recognized, the calibration information of this device will then be generated from the data base. Green for good to use, Red for out of calibration
- If the picture is blurry, device can not be recognized, an error message will be generated, so user can retake the picture



Mock-up --- 2) Select from Album


User can also take pictures from the album



Mock-up --- 3) Instruction

- The user should be able to use this app by intuition, but just in case, some instructions are aslo provided for the user:
 - Focus and put target in the center
 - Only picture one object at a time

Block Diagram



Expected Outcome

- With OneScan, quality assurance in North America can be improved in the long term:
 - the auditor can conduct more efficient/accurate auditing,
 - patients and car owners can have better quality products,
 - · factories will have less recalls
- North America's quality assurance system can then achieve the next level and will keep ahead in the Fourth Industrial Revolution (Industry4.0).





The End. Thank you!

